

Product datasheet

Anti-TRPM6 antibody ab47017

3 References 1 图像

概述

| | |
|-------|---|
| 产品名称 | Anti-TRPM6抗体 |
| 描述 | 豚鼠多克隆抗体to TRPM6 |
| 宿主 | Guinea pig |
| 经测试应用 | 适用于: IHC-Fr |
| 种属反应性 | 与反应: Rat 预测可用于: Mouse |
| 免疫原 | Synthetic peptide: GCERDKNRSSLERHTRL |

Run BLAST with

Run BLAST with

性能

| | |
|------|--|
| 形式 | Liquid |
| 存放说明 | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. |
| 存储溶液 | Preservative: 0.05% Sodium azide Constituents: PBS, Whole serum |
| 纯度 | Whole antiserum |
| 克隆 | 多克隆 |
| 同种型 | IgG |

应用

Our [Abpromise guarantee](#) covers the use of **ab47017** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


| 应用 | Ab评论 | 说明 |
|--------|------|-----------------|
| IHC-Fr | | 1/100 - 1/1000. |

靶标

相关性 TRPM6 is an essential ion channel and serine/threonine-protein kinase, and is crucial for magnesium homeostasis. TRPM6 also has an important role in epithelial magnesium transport and in the active magnesium absorption in the gut and kidney. The various isoforms of the type M6-kinase lack the ion channel region.

细胞定位 Multi-pass membrane protein.

图片



Immunohistochemical staining of frozen rat kidney tissue, using ab47017. Secondary staining was carried out using a Cy3 conjugated antibody.

Immunohistochemistry (Frozen sections) - Anti-TRPM6 antibody (ab47017)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors