

Product datasheet

Anti-UCK2 antibody ab60221

1 图像

概述

产品名称	Anti-UCK2抗体
描述	兔多克隆抗体to UCK2
宿主	Rabbit
经测试应用	适用于: ELISA, IHC-P
种属反应性	与反应: Human
免疫原	KLH conjugated synthetic peptide selected from the N-terminal region of human UCK2.
阳性对照	Human hepatocarcinoma tissue.

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.09% Sodium Azide Constituents: PBS
纯度	Protein G purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab60221** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
ELISA		1/1000.
IHC-P		1/50 - 1/100.

## 靶标


---

**相关性** UCK2 phosphorylates uridine and cytidine to uridine monophosphate and cytidine monophosphate. This is the first step in the production of the pyrimidine nucleoside triphosphates required for RNA and DNA synthesis. Does not phosphorylate deoxyribonucleosides or purine ribonucleosides.

**细胞定位** Cytosol

## 图片

---



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-UCK2 antibody (ab60221)

Ab60221 (1:50) staining human UCK2 in human hepatocarcinoma tissue by immunohistochemistry using formalin fixed, paraffin embedded tissue.

Formalin-fixed and paraffin-embedded human cancer tissue reacted with the primary antibody, which was peroxidase-conjugated to the secondary antibody, followed by AEC staining.

HC = hepatocarcinoma

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors