

## Product datasheet

# Anti-USP36 antibody ab80260

### 1 Abreviews

#### 概述

产品名称	Anti-USP36抗体
描述	兔多克隆抗体to USP36
宿主	Rabbit
特异性	Reacts specifically with human 122 kDa USP36 protein
经测试应用	适用于: IHC-Fr, IHC-P, WB
种属反应性	与反应: Human 预测可用于: Mouse 
免疫原	Synthetic peptide derived from internal domain of human USP36 protein

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
存储溶液	Preservative: None Constituents: Whole serum
纯度	Whole antiserum
克隆	多克隆
同种型	IgG

#### 应用

Our [Abpromise guarantee](#) covers the use of **ab80260** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
IHC-Fr		1/200 - 1/1000.
IHC-P		1/200 - 1/1000.

应用	Ab评论	说明
WB		1/500 - 1/5000. Detects a band of approximately 122 kDa.

## 靶标

功能	May be required for maintaining multiple types of adult stem cells. May function as a transcriptional repressor by continually deubiquitinating histone H2B at the promoters of genes critical for cellular differentiation, thereby preventing histone H3 'Lys-4' trimethylation (H3K4).
组织特异性	Broadly expressed.
序列相似性	Belongs to the peptidase C19 family.
细胞定位	Nucleus > nucleolus.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors