

Product datasheet

Anti-Ultraspiracle antibody ab106341

2 References

概述

产品名称	Anti-Ultraspiracle抗体
描述	兔多克隆抗体to Ultraspiracle
宿主	Rabbit
经测试应用	适用于: WB, ICC
种属反应性	与反应: Drosophila melanogaster
免疫原	Synthetic peptide derived from the N terminal domain of Ultraspiracle protein from Choristoneura fumiferana (Spruce budworm).

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C long term.
存储溶液	Constituent: Whole serum
纯度	Whole antiserum
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab106341** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/200 - 1/2000. Predicted molecular weight: 52 kDa.
ICC		Use at an assay dependent dilution.

靶标

功能 Receptor for ecdysone. May be an important modulator of insect metamorphosis. Plays an

important part in embryonic and post-embryonic development. Binds to ecdysone response elements (ECRES) such as in the promoter region of s15 chorion gene.

**序列相似性**

Belongs to the nuclear hormone receptor family. NR2 subfamily.  
Contains 1 nuclear receptor DNA-binding domain.

**结构域**

Composed of three domains: a modulating N-terminal domain, a DNA-binding domain and a C-terminal ligand-binding domain.

**细胞定位**

Nucleus.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

**Our Abpromise to you: Quality guaranteed and expert technical support**

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

**Terms and conditions**

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors