

Product datasheet

Anti-Vibrio parahaemolyticus K28 antibody ab79497

概述

产品名称	Anti-Vibrio parahaemolyticus K28抗体
描述	兔多克隆抗体to Vibrio parahaemolyticus K28
宿主	Rabbit
特异性	Specific to Vibrio parahaemolyticus capsular K antigen.
经测试应用	适用于: Agglutination
种属反应性	Vibrio parahaemolyticus K28
免疫原	Tissue/ cell preparation of Vibrio parahaemolyticus K28 (heat inactivated)
阳性对照	Antigenic suspension.
常规说明	Reagents should be allowed to stand at 15°C - 25°C for at least 30 minutes prior to use.

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.08% Sodium azide Constituent: Whole serum
纯度	Whole antiserum
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab79497** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
Agglutination		Use at an assay dependent concentration. See attached protocol.

靶标

相关性 *Vibrio parahaemolyticus* is a non-spore forming gram negative bacillus with monopolar flagella. This bacterium is frequently isolated from estuarine and marine environments and is the leading cause of seafood-associated gastroenteritis in the United States. Serological types are determined in combination with somatic antigens.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors