

Product datasheet

Anti-Vibrio parahaemolyticus K55 antibody ab78959

概述

产品名称	Anti-Vibrio parahaemolyticus K55抗体
描述	兔多克隆抗体to Vibrio parahaemolyticus K55
宿主	Rabbit
特异性	ab78959 is specific to Vibrio parahaemolyticus K55 antigen
经测试应用	适用于: Agglutination
种属反应性	Reacts with Vibrio parahaemolyticus K55
免疫原	Tissue/ cell preparation of Vibrio parahaemolyticus K55 (heat inactivated)
阳性对照	Antigenic suspension.
常规说明	Reagents should be allowed to stand at 15°C - 25°C for at least 30 minutes prior to use.

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.08% Sodium Azide Constituents: Whole serum
纯度	Whole antiserum
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab78959** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
Agglutination		Use at an assay dependent concentration. See attached protocol.

靶标

相关性 Vibrio parahaemolyticus, a member of the Vibrio genus, is a curved, rod shaped, Gram negative bacterium, found in brackish saltwater, which, when ingested, causes gastrointestinal illness in humans. Vibrio parahaemolyticus is oxidase positive, facultatively aerobic, non spore forming and is motile, with a single, polar flagellum. Most people become infected with Vibrio parahaemolyticus by eating raw or undercooked shellfish, particularly oysters. Less commonly, this organism can cause an infection in the skin when an open wound is exposed to warm seawater.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors