


Product datasheet

Anti-Wnt5b antibody ab115563

1 References 3 图像

概述

产品名称	Anti-Wnt5b抗体
描述	兔多克隆抗体to Wnt5b
宿主	Rabbit
经测试应用	适用于: IHC-P
种属反应性	与反应: Human 预测可用于: Mouse, Rat 
免疫原	Synthetic peptide conjugated to KLH, corresponding to internal sequence amino acids from Human Wnt5b (NP_116031.1).
阳性对照	Human Spleen, Human Thymus and Human Prostate tissues.

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at -80°C.
存储溶液	Preservative: 0.1% Sodium azide Constituent: 99% PBS
纯度	Protein A purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab115563** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
IHC-P		Use a concentration of 10 µg/ml. Formalin-fixed paraffin-embedded tissues requires pretreatment with Proteinase K.

脚标


---

<b>功能</b>	Ligand for members of the frizzled family of seven transmembrane receptors. Probable developmental protein. May be a signaling molecule which affects the development of discrete regions of tissues. Is likely to signal over only few cell diameters.
<b>序列相似性</b>	Belongs to the Wnt family.
<b>细胞定位</b>	Secreted > extracellular space > extracellular matrix.

---


## 图片

---




ab115563, at 10µg/ml, staining Wnt5b in Formalin-fixed, Paraffin-embedded Human Spleen tissue by Immunohistochemistry.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Wnt5b antibody (ab115563)



ab115563, at 10µg/ml, staining Wnt5b in Formalin-fixed, Paraffin-embedded Human Thymus tissue by Immunohistochemistry.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Wnt5b antibody (ab115563)



ab115563, at 10µg/ml, staining Wnt5b in Formalin-fixed, Paraffin-embedded Human Prostate tissue by Immunohistochemistry.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Wnt5b antibody (ab115563)

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors