


Product datasheet

Anti-Xenotropic MuLV-related virus antibody ab111014

1 图像

概述

产品名称	Anti-Xenotropic MuLV-related virus抗体
描述	兔多克隆抗体to Xenotropic MuLV-related virus
宿主	Rabbit
经测试应用	适用于: IHC-P
种属反应性	与反应: Human 预测可用于: Rat 
阳性对照	Prostate carcinoma with Xenotropic MuLV-related virus infection.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	pH: 7.60 Preservative: 0.1% Sodium azide Constituents: PBS, 1% BSA
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab111014** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab评论	说明
IHC-P		1/100. Boil tissue section in EDTA buffer for 10 min followed by cooling at RT for 20 min.

靶标

相关性

Xenotropic murine leukemia virus-related virus (XMRV) is a newly recognized human retrovirus first identified in prostate cancer tissues from men with a deficiency in the innate immunity gene RNASEL. At present, studies have detected XMRV DNA and protein in prostate cancer at widely different rates and in patients with chronic fatigue syndrome. The XMRV protein is synthesized in the Golgi by host cells.

图片



ab111014 at 1/100 dilution staining
Xenotropic MuLV-related virus in Formalin-fixed, paraffin-embedded Human prostate adenocarcinoma by Immunohistochemistry.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Xenotropic MuLV-related virus antibody (ab111014)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors