


## Product datasheet

# Anti-ZDHHC20 antibody ab110478

1 References 1 图像

### 概述

产品名称	Anti-ZDHHC20抗体
描述	兔多克隆抗体to ZDHHC20
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Human 预测可用于: Mouse, Rat 
免疫原	Synthesized peptide derived from an internal sequence of Human ZDHHC20.
阳性对照	Extracts from COLO, MCF7 and Jurkat cells

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at -20°C. Stable for 12 months at -20°C.
存储溶液	pH: 7.40 Preservative: 0.02% Sodium azide Constituents: 99% PBS, 0.88% Sodium chloride, 50% Glycerol Note: PBS is without Mg <sup>2+</sup> and Ca <sup>2+</sup> .
纯度	Immunogen affinity purified
纯化说明	ab110478 was affinity purified from rabbit antiserum by affinity chromatography using epitope-specific immunogen.
克隆	多克隆
同种型	IgG

### 应用

Our [Abpromise guarantee](#) covers the use of **ab110478** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/500 - 1/1000. Predicted molecular weight: 42 kDa.

## 靶标

序列相似性	Belongs to the DHHC palmitoyltransferase family. Contains 1 DHHC-type zinc finger.
结构域	The DHHC domain is required for palmitoyltransferase activity.
翻译后修饰	Phosphorylated upon DNA damage, probably by ATM or ATR.
细胞定位	Membrane.

## 图片



**All lanes :** Anti-ZDHHC20 antibody (ab110478) at 1/500 dilution

**Lane 1 :** Extracts from COLO cells

**Lane 2 :** Extracts from MCF7 cells

**Lane 3 :** Extracts from Jurkat cells

**Lane 4 :** Extracts from Jurkat cells with immunized peptide at 10 µg

Lysates/proteins at 30 µg per lane.

**Predicted band size:** 42 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors