


Product datasheet

Anti-ZNF346 antibody ab96198

2 图像

概述

产品名称	Anti-ZNF346抗体
描述	兔多克隆抗体to ZNF346
宿主	Rabbit
经测试应用	适用于: WB, IHC-P
种属反应性	与反应: Human 预测可用于: Mouse, Cow 
免疫原	Recombinant fragment, corresponding to a region within amino acids 1-274 of Human ZNF346 (NP_036411).
阳性对照	HeLa and HepG2 whole cell lysates; A549 xenograft; 293T, A431, H1299, Molt-4 and Raji cell lysates

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: 0.01% Thimerosal (merthiolate) Constituents: 10% Glycerol, 0.1M Tris, 0.1M Glycine, pH 7.0
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab96198** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/500 - 1/3000. Predicted molecular weight: 33 kDa.

应用	Ab评论	说明
IHC-P		1/100 - 1/500.

## 靶标

## 图片



**All lanes :** Anti-ZNF346 antibody (ab96198)  
at 1/5000 dilution


**Lane 1 :** HeLa whole cell lysate

**Lane 2 :** HepG2 whole cell lysate

Lysates/proteins at 30 µg per lane.

**Predicted band size:** 33 kDa

12% SDS PAGE



Immunohistochemical analysis of ZNF346 in paraffin-embedded A549 xenograft, using ab96198 at 1/100 dilution.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery

- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors