

Product datasheet

Anti-ZNF75 antibody [3424C1a] ab76563

1 图像

概述

产品名称	Anti-ZNF75抗体[3424C1a]
描述	小鼠单克隆抗体[3424C1a] to ZNF75
宿主	Mouse
经测试应用	适用于: WB, Dot blot
种属反应性	与反应: Human
免疫原	Recombinant fragment comprising of 50-200 amino acids from the internal region of Human ZNF75 (NM_007131).
阳性对照	Recombinant fragment of human ZNF75.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: 0.05% Sodium Azide Constituents: 1% BSA, PBS, (8.0mM Sodium phosphate, 3.0mM Potassium chloride, 140mM Sodium chloride, 1.5mM Potassium phosphate), pH 7.4
纯度	Protein G purified
纯化说明	ab76563 was purified using protein G column chromatography from culture supernatant of hybridoma cultured in a medium containing bovine IgG-depleted (approximately 95%) fetal bovine serum and filtered through a 0.22µm membrane.
克隆	单克隆
克隆编号	3424C1a
同种型	IgG1

应用

Our [Abpromise guarantee](#) covers the use of **ab76563** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab 评论	说明
WB		Use a concentration of 0.2 - 2 µg/ml. Predicted molecular weight: 59 kDa. A band of 32 kDa was detected when the immunizing recombinant protein fragment was used.
Dot blot		Use at an assay dependent dilution.

靶标

功能	May be involved in transcriptional regulation.
序列相似性	Belongs to the krueppel C2H2-type zinc-finger protein family. Contains 5 C2H2-type zinc fingers. Contains 1 KRAB domain. Contains 1 SCAN box domain.
细胞定位	Nucleus.

图片



Anti-ZNF75 antibody [3424C1a] (ab76563) at 2 µg/ml + immunizing recombinant protein

Predicted band size: 59 kDa

Observed band size: 32 kDa

The molecular weight of the band on the western blot does not correspond to the molecular weight of the natural protein because only a fragment of the protein was used.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors