

## Product datasheet

# Anti-Marburg Virus antibody [FM44] ab10007

### 概述

产品名称	Anti-Marburg Virus抗体[FM44]
描述	小鼠单克隆抗体[FM44] to Marburg Virus
宿主	Mouse
特异性	<a href="#">ab2055</a> , 2058, 2059, 10003, 10004, 10005, 10007 and 10009 recognise different epitopes which have not been fully determined.
经测试应用	<b>适用于:</b> ELISA, WB
免疫原	Marburg virus.
常规说明	Concentration varies from lot to lot and can be provided on request.

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: 0.1% Sodium Azide Constituents: PBS, pH 7.2
纯度	Protein G purified
纯化说明	Purity tested by elctrophoresis.
克隆	单克隆
克隆编号	FM44
同种型	IgG1

### 应用

Our [Abpromise guarantee](#) covers the use of **ab10007** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
ELISA		Use at an assay dependent dilution.
WB		Use at an assay dependent dilution. Detects a band of approximately 30 kDa.

## 靶标

---

**相关性** Marburg virus (MBG) is a representative of Filoviridae family of RNA containing viruses. MBG is an exceptionally dangerous pathogen, which induces a severe contagious and highly lethal (53-88%) febrile disease with hemorrhagic syndrome. The name “thread viruses” is based on their morphology. MBG made its appearance in 1967 in the form of frightening nosocomial outbreak, initially among polio vaccine production workers in Germany in contact with Ugandan green monkeys and their kidney tissues. The patho-and immuno-genesis of MBG fever have so far been little studied, no specific agents or treatment methods have been developed.

**细胞定位** Viral

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors