

## Product datasheet

# Anti-Rotavirus group specific antigen antibody [3C10] ab10012

### 概述

产品名称	Anti-Rotavirus group specific antigen抗体[3C10]
描述	小鼠单克隆抗体[3C10] to Rotavirus group specific antigen
宿主	Mouse
特异性	Rotavirus P42 inner-capsid antigen specific. Cross-reacts with monkey rotavirus (SA-11), porcine rotavirus (PP) and human rotaviruses (field strains).
经测试应用	适用于: WB
免疫原	Full length native protein (purified bovine rotavirus, cell culture adapted strain MR).
常规说明	Concentration varies from lot to lot and can be provided on request.

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: 0.1% Sodium Azide Constituents: PBS, pH 7.4
纯度	Protein A purified
纯化说明	Purity tested by elcetrophoresis.
克隆	单克隆
克隆编号	3C10
骨髓瘤	Sp2/0
同种型	IgG2a

### 应用

Our [Abpromise guarantee](#) covers the use of **ab10012** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use at an assay dependent dilution.

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors