

Product datasheet

Anti-Nkx2.4 antibody ab105826

1 图像

概述

产品名称	Anti-Nkx2.4抗体
描述	兔多克隆抗体to Nkx2.4
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Mouse 预测可用于: Rat, Dog, Human, Pig 
免疫原	Synthetic peptide corresponding to a region within the C terminal amino acids 301-350 (QQPQAPTAP ELEELSPSP ALHGPGGGLA ALDAATGDY GGVLGANLLY) of Mouse Nkx2.4 (NP_075993). <a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a>
阳性对照	Testis lysate

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: 0.09% Sodium azide Constituents: 2% Sucrose, PBS
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab105826** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab评论	说明

应用	Ab评论	说明
WB		Use a concentration of 0.5 µg/ml. Predicted molecular weight: 36 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

## 靶标

功能	Probable transcription factor.
序列相似性	Belongs to the NK-2 homeobox family. Contains 1 homeobox DNA-binding domain.
细胞定位	Nucleus.

## 图片



Anti-Nkx2.4 antibody (ab105826) at 0.5 µg/ml + Mouse testis lysate at 10 µg

**Predicted band size:** 36 kDa

Gel concentration: 12%

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors