

Product datasheet

Anti-Proteasome 19S S5A antibody [21] ab106186

1 图像

概述

产品名称	Anti-Proteasome 19S S5A抗体[21]
描述	小鼠单克隆抗体[21] to Proteasome 19S S5A
宿主	Mouse
经测试应用	适用于: WB, ELISA
种属反应性	与反应: Human
免疫原	Recombinant protein corresponding to full length Human Proteasome 19S S5A expressed in E. coli.
阳性对照	Jurkat whole cell lysate

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
存储溶液	Preservative: 0.02% Proclin Constituents: Tris HCl, Sodium citrate, pH 7.0
纯度	Protein G purified
克隆	单克隆
克隆编号	21
同种型	IgG1

应用

Our [Abpromise guarantee](#) covers the use of **ab106186** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/1000. Predicted molecular weight: 41 kDa.
ELISA		Use at an assay dependent dilution.

靶标

功能	Binds and presumably selects ubiquitin-conjugates for destruction. Displays selectivity for longer polyubiquitin chains. Modulates intestinal fluid secretion.
序列相似性	Belongs to the proteasome subunit S5A family. Contains 2 UIM (ubiquitin-interacting motif) repeats. Contains 1 VWFA domain.
结构域	The 2 UIM motifs are involved in the binding to a multi-ubiquitin chain in a cooperative way.
翻译后修饰	Phosphorylated upon DNA damage, probably by ATM or ATR.

图片



Anti-Proteasome 19S S5A antibody [21]
(ab106186) at 1/1000 dilution + Jurkat whole
cell lysate

Predicted band size: 41 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors