

Product datasheet

Anti-DDDDK tag antibody [29E4.G7] - (Equivalent to FLAG antibodies from Sigma) (FITC) ab106221

概述

产品名称	Anti-DDDDK tag抗体[29E4.G7] - (Equivalent to FLAG抗体from Sigma) (FITC)
描述	小鼠单克隆抗体[29E4.G7] to DDDDK tag - (Equivalent to FLAG抗体from Sigma) (FITC)
宿主	Mouse
偶联物	FITC. Ex: 493nm, Em: 528nm
特异性	ab106221 recognizes the DYKDDDDK epitope tag (Asp-Tyr-Lys-Asp-Asp-Asp-Asp-Lys) fused to either the amino- or carboxy- termini of targeted proteins in transfected or transformed cells.
免疫原	Synthetic peptide corresponding to DDDDK tag conjugated to Keyhole Limpet Haemocyanin (KLH).

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 7.20 Preservative: 0.01% Sodium azide Constituents: 0.42% Potassium phosphate, 0.88% Sodium chloride, 1% BSA
纯度	Ascites
克隆	单克隆
克隆编号	29E4.G7
同种型	IgG

靶标

相关性	This is a useful tool for the localisation and characterisation of DDDDK tagged proteins.
------------	---

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors