


Product datasheet

Anti-SOX22 antibody ab106339

1 图像

概述

产品名称	Anti-SOX22抗体
描述	兔多克隆抗体to SOX22
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Human 预测可用于: Mouse, Rat, Horse, Guinea pig, Dog, Pig 
免疫原	Synthetic peptide corresponding to a region within C terminal amino acids 266-315 (DCSALDRDPD LQPPSGTSHF EFPDYCTPEV TEMIAGDWRP SSIADLVFTY) of Human SOX22 (NP_008874). Run BLAST with ExPASy Run BLAST with NCBI
阳性对照	Human fetal brain lysate.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: 0.09% Sodium azide Constituents: 2% Sucrose, PBS
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab106339** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab评论	说明

应用	Ab评论	说明
WB		Use a concentration of 0.25 µg/ml. Predicted molecular weight: 34 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

靶标

功能	Binds to the sequence 5'-AACAAAT-3'.
组织特异性	Expressed most abundantly in the CNS. Also expressed in fetal brain and kidney and adult heart, pancreas, testis and ovary. Other tissues were only weakly positive.
序列相似性	Contains 1 HMG box DNA-binding domain.
细胞定位	Nucleus.

图片



Anti-SOX22 antibody (ab106339) at 0.25 µg/ml + Human fetal brain lysate at 10 µg

Predicted band size: 34 kDa

Gel concentration: 12%.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors