

Product datasheet

Anti-GPR83 antibody ab108432

1 图像

概述

|       |  |
|-------|--|
| 产品名称  | Anti-GPR83抗体   |
| 描述    | 兔多克隆抗体to GPR83   |
| 宿主    | Rabbit   |
| 经测试应用 | 适用于: WB  |
| 种属反应性 | 与反应: Human   |
| 免疫原   | Synthetic peptide corresponding to a region within C terminal amino acids 370-419 (QEDRPPSPVP SFRVAWTEKN DGQRAPLANN LLPTSQLQSG KTDLSSVEPI) of Human GPCR GPR72 (NP_057624).<br><a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a> |
| 阳性对照  | 721_B cell lysate.   |

性能

|      |  |
|------|--|
| 形式   | Liquid   |
| 存放说明 | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles. |
| 存储溶液 | Constituents: PBS, 2% Sucrose  |
| 纯度   | Immunogen affinity purified  |
| 克隆   | 多克隆  |
| 同种型  | IgG  |

应用

Our [Abpromise guarantee](#) covers the use of **ab108432** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| 应用 | Ab评论 | 说明   |
|----|------|--|
| WB |      | Use a concentration of 1 µg/ml. Predicted molecular weight: 48 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T. |

## 靶标

|       |  |
|-------|--|
| 功能    | Orphan receptor. Could be a neuropeptide Y receptor. |
| 组织特异性 | Brain specific.                                      |
| 序列相似性 | Belongs to the G-protein coupled receptor 1 family.  |
| 细胞定位  | Cell membrane.                                       |

## 图片



Anti-GPR83 antibody (ab108432) at 1  $\mu$ g/ml +  
721\_B cell lysate at 10  $\mu$ g

**Predicted band size:** 48 kDa

Gel concentration: 12%.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors