

Product datasheet

Anti-B4GALT5 antibody ab110720

1 图像

概述

| | |
|-------|---|
| 产品名称 | Anti-B4GALT5抗体 |
| 描述 | 兔多克隆抗体to B4GALT5 |
| 宿主 | Rabbit |
| 经测试应用 | 适用于: WB |
| 种属反应性 | 与反应: Mouse 预测可用于: Rat, Rabbit, Horse, Chicken, Guinea pig, Cow, Cat, Dog, Pig, Zebrafish |
| 免疫原 | Synthetic peptide corresponding to a region within C-terminal amino acids 311-360 (GYSVSRPEGD TGKYKSIPHH HRGEVQFLGR YALLRKSKER QGLDGLNNLN) of Mouse B4GALT5 (NP_062809). Run BLAST with ExPASy Run BLAST with NCBI |
| 阳性对照 | Mouse Kidney Lysate. |

性能

| | |
|------|--|
| 形式 | Liquid |
| 存放说明 | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles. |
| 存储溶液 | Preservative: 0.09% Sodium azide Constituents: 2% Sucrose, PBS |
| 纯度 | Immunogen affinity purified |
| 克隆 | 多克隆 |
| 同种型 | IgG |

应用

Our [Abpromise guarantee](#) covers the use of **ab110720** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


| 应用 | Ab评论 | 说明 |
|----|------|----|
| | | |

| 应用 | Ab评论 | 说明 |
|----|------|--|
| WB | | Use a concentration of 1 µg/ml. Predicted molecular weight: 45 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T. |

靶标

| | |
|-------|--|
| 功能 | Responsible for the synthesis of complex-type N-linked oligosaccharides in many glycoproteins as well as the carbohydrate moieties of glycolipids. |
| 组织特异性 | Ubiquitously expressed. |
| 通路 | Protein modification; protein glycosylation. |
| 序列相似性 | Belongs to the glycosyltransferase 7 family. |
| 细胞定位 | Golgi apparatus > Golgi stack membrane. Trans cisternae of Golgi stack. |

图片



Anti-B4GALT5 antibody (ab110720) at 1 µg/ml + Mouse Kidney
Lysate at 10 µg

Predicted band size: 45 kDa

12% SDS PAGE

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise,

please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors