

## Product datasheet

# Anti-MEKK3 antibody ab111259

### 1 图像

#### 概述

产品名称	Anti-MEKK3抗体
描述	兔多克隆抗体to MEKK3
宿主	Rabbit
经测试应用	适用于: WB, IHC-P
种属反应性	与反应: Mouse 预测可用于: Rat, Sheep, Rabbit, Guinea pig, Hamster, Dog, Human, Pig, Xenopus laevis, Drosophila melanogaster, Monkey
免疫原	Synthetic peptide derived from sequence near the N-terminus of Mouse MEKK3.
阳性对照	Human Uterus tissue

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Constituent: Whole serum
纯度	Whole antiserum
克隆	多克隆
同种型	IgG

#### 应用

Our [Abpromise guarantee](#) covers the use of **ab111259** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/5000. Predicted molecular weight: 71 kDa.
IHC-P		1/100. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.

## 靶标


---

<b>功能</b>	Component of a protein kinase signal transduction cascade. Mediates activation of the NF-kappa-B, AP1 and DDIT3 transcriptional regulators.
<b>序列相似性</b>	Belongs to the protein kinase superfamily. STE Ser/Thr protein kinase family. MAP kinase kinase subfamily. Contains 1 OPR domain. Contains 1 protein kinase domain.

---

## 图片

---



Formalin-Fixed, Paraffin-Embedded Human Uterus tissue stained using ab111259 at a dilution of 1/100.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-MEKK3 antibody (ab111259)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors