

## Product datasheet

# Anti-FAM53A antibody ab111497

### 1 图像

#### 概述

产品名称	Anti-FAM53A抗体
描述	兔多克隆抗体to FAM53A
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Synthetic peptide, corresponding to a region within C terminal amino acids 341-390 (RGLRTSPVHP NLWASRESVT SDGSRRSSGD PRDGDSVGEE GVFPRARWEL) of Human FAM53A (NP_001013644). <a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a>
阳性对照	HeLa cells lysate

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Constituents: PBS, 2% Sucrose
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

#### 应用

Our [Abpromise guarantee](#) covers the use of **ab111497** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab评论	说明
----	------	----

WB		Use a concentration of 0.125 µg/ml. Predicted molecular weight: 43 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.
----	--	--

## 靶标

功能	May play an important role in neural development; the dorsomedial roof of the third ventricle.
序列相似性	Belongs to the FAM53 family.
细胞定位	Nucleus. Subnuclear distribution.

## 图片



Anti-FAM53A antibody (ab111497) at 0.125  $\mu$ g/ml + HeLa cell lysate at 10  $\mu$ g

**Predicted band size:** 43 kDa

Gel concentration: 12%

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors