

Product datasheet

Anti-T Cell Receptor antibody [KJ1-26] (FITC) ab114058

概述

| | |
|-------|---|
| 产品名称 | Anti-T Cell Receptor抗体[KJ1-26] (FITC) |
| 描述 | 小鼠单克隆抗体[KJ1-26] to T Cell Receptor (FITC) |
| 宿主 | Mouse |
| 偶联物 | FITC. Ex: 493nm, Em: 528nm |
| 特异性 | ab114058 binds to the T-cell receptor (TCR) expressed on DO11.10 and DO11.10.24 Murine T cell hybridomas, but not on any other T cell hybridomas. |
| 经测试应用 | 适用于: Flow Cyt |
| 种属反应性 | 与反应: Mouse |
| 免疫原 | T cell hybridoma, D011.10 |
| 阳性对照 | Balb/c T Cell Receptor Positive Thymocytes |

性能

| | |
|------|---|
| 形式 | Liquid |
| 存放说明 | Shipped at 4°C. Store at +4°C. |
| 存储溶液 | Preservative: 0.02% Sodium azide Constituents: 0.5% BSA, 99% PBS |
| 克隆 | 单克隆 |
| 克隆编号 | KJ1-26 |
| 同种型 | IgG2a |

应用

Our [Abpromise guarantee](#) covers the use of **ab114058** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| 应用 | Ab评论 | 说明 |
|----------|------|---|
| Flow Cyt | | Use 0.5µg for 10 ⁶ cells. ab91362 - Mouse monoclonal IgG2a, is suitable for use as an isotype control with this antibody. |

靶标

相关性 Two distinct types of T-cell antigen receptors have been identified: the alpha/beta heterodimer found on functional helper and cytotoxic T cells, and the gamma/delta heterodimer. The latter is first detected approximately 2 days before the appearance of cell-surface alpha/beta heterodimer during T-cell ontogeny.

细胞定位 Cell Membrane

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors