


Product datasheet

Anti-DNAJC11 antibody ab118744

1 图像

概述

产品名称	Anti-DNAJC11抗体
描述	兔多克隆抗体to DNAJC11
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Human 预测可用于: Mouse, Rat, Rabbit, Horse, Chicken, Guinea pig, Cow, Cat, Dog, Saccharomyces cerevisiae, Zebrafish 
免疫原	Synthetic peptide, corresponding to a region within C terminal amino acids 407-456 (IIKPYLRAQK EKELEKQRES AATDVLQKKQ EAESAVRLMQ ESVRRIIIEAE) of Human DNAJC11 (NP_060668.2). Run BLAST with ExPASy Run BLAST with NCBI
阳性对照	Human fetal heart lysate.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	Constituents: 97% PBS, 2% Sucrose
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab118744** in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明


应用	Ab评论	说明
----	------	----

WB Use a concentration of 1 µg/ml. Predicted molecular weight: 63 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

靶标

序列相似性 Belongs to the DNAJC11 family.
Contains 1 J domain.

图片



Anti-DNAJC11 antibody (ab118744) at 1 µg/ml + Human fetal heart lysate at 10 µg

Predicted band size: 63 kDa

Gel concentration: 12%

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors