

## Product datasheet

# Anti-FAM119A antibody [2D6] ab119186

### 1 图像

#### 概述

产品名称	Anti-FAM119A抗体[2D6]
描述	小鼠单克隆抗体[2D6] to FAM119A
宿主	Mouse
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Recombinant full length Human FAM119A (NP_660323) produced in HEK293T cell.
阳性对照	HEK293T cells transfected with the pCMV6-ENTRY FAM119A

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	pH: 7.30 Preservative: 0.02% Sodium azide Constituents: 48% PBS, 1% BSA, 50% Glycerol
纯度	IgG fraction
纯化说明	Purified from mouse ascites fluids by affinity chromatography
克隆	单克隆
克隆编号	2D6
同种型	IgG1

#### 应用

Our [Abpromise guarantee](#) covers the use of **ab119186** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/2000. Predicted molecular weight: 25 kDa.

## 靶标

序列相似性	Belongs to the FAM119 family.
细胞定位	Membrane.

## 图片



**All lanes** : Anti-FAM119A antibody [2D6]  
(ab119186) at 1/2000 dilution

**Lane 1** : HEK293T cells transfected with the  
pCMV6-ENTRY control

**Lane 2** : HEK293T cells transfected with the  
pCMV6-ENTRY FAM119A

Lysates/proteins at 5 µg per lane.

**Predicted band size:** 25 kDa

HEK293T cell lysates were generated from  
transient transfection of the cDNA clone  
(RC202889)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors