

## Product datasheet

# Anti-Influenza B Virus antibody [8.F.183] ab119967

### 概述

产品名称	Anti-Influenza B Virus抗体[8.F.183]
描述	小鼠单克隆抗体[8.F.183] to Influenza B Virus
宿主	Mouse
特异性	ab119967 recognizes Influenza B nucleoprotein.
经测试应用	适用于: WB, ELISA, Flow Cyt, IHC-P
种属反应性	与反应: Other species
免疫原	Influenza B Virus, Singapore strain.

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Add glycerol to a final volume of 50% for extra stability and aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
存储溶液	pH: 7.30 Preservative: 0.05% Sodium azide Constituent: 99% PBS
纯度	Protein A purified
克隆	单克隆
克隆编号	8.F.183
同种型	IgG2

### 应用

Our [Abpromise guarantee](#) covers the use of **ab119967** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use at an assay dependent concentration.
ELISA		Use at an assay dependent concentration.

应用	Ab评论	说明
Flow Cyt		Use at an assay dependent concentration. <a href="#">ab18414</a> - Mouse monoclonal IgG2a, is suitable for use as an isotype control with this antibody.
IHC-P		Use at an assay dependent concentration.

## 靶标

相关性	Influenza A and B are the two types of influenza viruses that cause epidemic human disease. Influenza B viruses currently circulating are divided into two antigenically and genetically distinct groups: the Victoria and Yamagata lineages.
细胞定位	Nuclear

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors