

Anti-Intra Acrosomal Protein antibody [Ds-2] ab14300

概述

产品名称	Anti-Intra Acrosomal蛋白抗体[Ds-2]
描述	小鼠单克隆抗体[Ds-2] to Intra Acrosomal蛋白
宿主	Mouse
经测试应用	适用于: WB, ICC
种属反应性	与反应: Dog
免疫原	Tissue/ cell preparation (Dog). Freshly ejaculated dog sperms were washed in PBS and extracted in 3% acetic acid, 10% glycerol, 30 mM benzaminidine. The acid extract was dialyzed against 0.2% acetic acid, lyophilized and subsequently used for immunization.

常规说明

The antibody reacts with intra-acrosomal proteins (92 and 42 kDa) of canine sperms.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: 15 mM sodium azide Constituents: Phosphate buffered saline, pH 7.4
纯度	>95% by SDS-PAGE
Primary antibody说明	The antibody reacts with intra-acrosomal proteins (92 and 42 kDa) of canine sperms.
克隆	单克隆
克隆编号	Ds-2
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab14300** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use at an assay dependent concentration. Predicted molecular weight: 32 kDa.
ICC		Use at an assay dependent concentration. Membrane permeabilization (acetone) essential.

靶标

组织特异性	Testis.
细胞定位	Cytoplasmic vesicle > secretory vesicle > acrosome. Nascent acrosomal vesicle of Golgi phase spermatids.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors