

Product datasheet

Anti-Transmissible Gastroenteritis Virus antibody [8004] ab20301

概述

产品名称	Anti-Transmissible Gastroenteritis Virus抗体[8004]
描述	小鼠单克隆抗体[8004] to Transmissible Gastroenteritis Virus
宿主	Mouse
特异性	Recognizes the peplomer protein of TGEV.
经测试应用	适用于: Neutralising
种属反应性	Does not react with the following coronaviruses: FIPV, FeCV, and CCV. Not yet tested in other species.
免疫原	Tissue/ cell preparation (TGEV-infected ST cells).

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: None Constituents: PBS, pH 7.4
纯度	Protein G purified
克隆	单克隆
克隆编号	8004
骨髓瘤	unknown
同种型	IgG2a
轻链类型	kappa

应用

Our [Abpromise guarantee](#) covers the use of **ab20301** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
Neutralising		

应用说明

Neut: Use at an assay dependent dilution. May be used in immunoassays to identify TGEV as well as TGEV neutralization assays.

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

靶标

相关性

Transmissible gastroenteritis (TGE) is an acute highly contagious disease of pigs caused by virus from the Coronaviridae family, Coronavirus genus. The disease is characterised by profuse diarrhoea and vomiting. There is high morbidity and mortality in piglets.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors