

Product datasheet

Anti-Yersinia pestis v antigen antibody [Va48] ab20983

概述

产品名称	Anti-Yersinia pestis v antigen抗体[Va48]
描述	小鼠单克隆抗体[Va48] to Yersinia pestis v antigen
宿主	Mouse
特异性	Yersinia pestis, V antigen
经测试应用	适用于: ELISA, WB
种属反应性	Reacts with Yersinia pestis.
免疫原	Recombinant full length protein

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: 0.1% Sodium Azide Constituents: PBS, pH 7.4
纯度	Protein G purified
克隆	单克隆
克隆编号	Va48
同种型	IgG1

应用

Our [Abpromise guarantee](#) covers the use of **ab20983** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
ELISA		Use at an assay dependent dilution.
WB		Use at an assay dependent dilution. Predicted molecular weight: 37 kDa.

图标

相关性	Yerisinia pestis has a number of virulence factors that enable it to survive in humans by facilitating use of host nutrients, causing damage to host cells, and subverting phagocytosis and other host defense mechanisms. Yerisinia pestis is possibly involved in calcium regulation of YOP expression, which includes the export process.
细胞定位	Secreted

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors