

## Product datasheet

# Anti-CD62L antibody [LT-TD180] ab21907

### 3 References

#### 概述

产品名称	Anti-CD62L抗体[LT-TD180]
描述	小鼠单克隆抗体[LT-TD180] to CD62L
宿主	Mouse
特异性	Recognises CD62L (L selectin)
经测试应用	适用于: Flow Cyt
种属反应性	与反应: Human
免疫原	Tissue/ cell preparation (Human): Peripheral blood leukocytes.

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 15mM Sodium Azide Constituents: PBS, pH 7.4
纯度	>95% by SDS-PAGE
纯化说明	Purified from ascites by DEAE chromatography and precipitation methods.
克隆	单克隆
克隆编号	LT-TD180
同种型	IgG1

#### 应用

Our [Abpromise guarantee](#) covers the use of **ab21907** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
----	------	----

应用	Ab评论	说明
----	------	----

Flow Cyt

Use at an assay dependent concentration.  
[ab170190](#) - Mouse monoclonal IgG1, is suitable for use as an isotype control with this antibody.

## 靶标

<b>功能</b>	Cell surface adhesion protein. Mediates the adherence of lymphocytes to endothelial cells of high endothelial venules in peripheral lymph nodes. Promotes initial tethering and rolling of leukocytes in endothelia.
<b>组织特异性</b>	Expressed in B cell lines and T lymphocytes.
<b>序列相似性</b>	Belongs to the selectin/LECAM family. Contains 1 C-type lectin domain. Contains 1 EGF-like domain. Contains 2 Sushi (CCP/SCR) domains.
<b>细胞定位</b>	Membrane.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors