

Product datasheet

Anti-Interferon alpha 1 antibody [4E-A1] ab24311

1 References

概述

产品名称	Anti-Interferon alpha 1抗体[4E-A1]
描述	大鼠单克隆抗体[4E-A1] to Interferon alpha 1
宿主	Rat
特异性	This antibody reacts with mouse interferon alpha (alpha 1, alpha 2 and alpha A), but not mouse interferon beta or gamma.
经测试应用	适用于: Neutralising
种属反应性	与反应: Mouse 不与反应: Human
免疫原	Purified mouse interferon alpha.
常规说明	Abcam is committed to meeting high standards of ethical manufacturing and as such, we will be discontinuing this product, which has been generated by the ascites method, within the next year. We are sorry for any inconvenience this may cause. If you would like help finding an alternative product, please do not hesitate to contact our scientific support team.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None. Constituents: PBS, pH 7.4.
纯度	Protein G purified
克隆	单克隆
克隆编号	4E-A1
同种型	IgG1

应用

Our [Abpromise guarantee](#) covers the use of **ab24311** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
Neutralising		Use at an assay dependent dilution.

靶标

功能	Produced by macrophages, IFN-alpha have antiviral activities. Interferon stimulates the production of two enzymes: a protein kinase and an oligoadenylate synthetase.
序列相似性	Belongs to the alpha/beta interferon family.
细胞定位	Secreted.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors