

## Product datasheet

# Anti-Ly6g antibody [RB6-8C5] (Phycoerythrin) ab24884

### 11 References

#### 概述

产品名称	Anti-Ly6g抗体[RB6-8C5] (Phycoerythrin)
描述	大鼠单克隆抗体[RB6-8C5] to Ly6g (Phycoerythrin)
宿主	Rat
偶联物	Phycoerythrin. Ex: 488nm, Em: 575nm
经测试应用	适用于: Flow Cyt
种属反应性	与反应: Mouse
免疫原	The details of the immunogen for this antibody are not available.

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at +4°C. Do Not Freeze.
存储溶液	pH: 7.20 Preservative: 0.09% Sodium azide Constituents: 0.1% Gelatin, 0.87% Sodium chloride, 0.142% Monobasic dihydrogen sodium phosphate
纯度	IgG fraction
纯化说明	Purified from TCS
克隆	单克隆
克隆编号	RB6-8C5
同种型	IgG2b
轻链类型	kappa

#### 应用

Our [Abpromise guarantee](#) covers the use of **ab24884** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
Flow Cyt		Use at an assay dependent concentration.

## 靶标

<b>相关性</b>	Ly6G is a GPI-anchored protein, that is also known as the myeloid differentiation antigen Gr1. The antigen is transiently expressed on monocytes in the bone marrow. The level of antigen expression in the bone marrow directly correlates with granulocyte differentiation and maturation. Ly6G is expressed predominantly on neutrophils, also in a subset of eosinophils, differentiating pre-monocytes, and plasmacytoid dendritic cells.
<b>细胞定位</b>	Attached to the membrane by a GPI-anchor.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors