

Product datasheet

Anti-MHC Class I H2 Dk antibody [15-5-5.3] (FITC) ab25142

概述

产品名称	Anti-MHC Class I H2 Dk抗体[15-5-5.3] (FITC)
描述	小鼠单克隆抗体[15-5-5.3] to MHC Class I H2 Dk (FITC)
宿主	Mouse
偶联物	FITC. Ex: 493nm, Em: 528nm
特异性	MAb 15-5-5.3 reacts with the H2D ^k class I alloantigen. It cross-reacts with H2K ^d and with cells from mice with the H2 ^f haplotype. The antibody does not react with other (e.g., a, b, p, q, r, s) haplotypes.
经测试应用	适用于: Flow Cyt
种属反应性	与反应: Mouse
免疫原	Tissue/ cell preparation of splenocytes (Mouse)

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: 0.1% Sodium Azide Constituents: PBS
纯度	Protein A purified
克隆	单克隆
克隆编号	15-5-5.3
同种型	IgG2a
轻链类型	kappa

应用

Our [Abpromise guarantee](#) covers the use of **ab25142** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
----	------	----

应用	Ab评论	说明
----	------	----

Flow Cyt

Use 1µg for 10⁶ cells.

[ab18455](#) - Mouse monoclonal IgG2a, is suitable for use as an isotype control with this antibody.

靶标

相关性 MHC Class I molecules play a central role in the immune system by presenting peptides derived from the endoplasmic reticulum lumen. MHC class I antigens are heterodimers consisting of one alpha chain (44kDa) with beta 2 microglobulin (11.5 kDa).

细胞定位 Cell Membrane

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors