

## Product datasheet

# Anti-RLA DQ antibody [2C4] (FITC) ab25685

### 2 References

#### 概述

产品名称	Anti-RLA DQ抗体[2C4] (FITC)
描述	小鼠单克隆抗体[2C4] to RLA DQ (FITC)
宿主	Mouse
偶联物	FITC. Ex: 493nm, Em: 528nm
特异性	ab25685 reacts with RLA-DQ-transfected cells, but not with RLA-DR-transfected cells.
经测试应用	<b>适用于:</b> Flow Cyt, IHC-Fr, IP
种属反应性	<b>与反应:</b> Rabbit
免疫原	Tissue/ cell preparation (Rabbit) (spleen cells)

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: 0.1% Sodium Azide Constituents: PBS
纯度	IgG fraction
克隆	单克隆
克隆编号	2C4
同种型	IgG2a
轻链类型	kappa

#### 应用

Our [Abpromise guarantee](#) covers the use of **ab25685** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
----	------	----

应用	Ab评论	说明
Flow Cyt		Use 1µg for 10 <sup>6</sup> cells. <a href="#">ab18455</a> - Mouse monoclonal IgG2a, is suitable for use as an isotype control with this antibody.
IHC-Fr		Use at an assay dependent concentration.
IP		Use at an assay dependent concentration.

## 靶标

**相关性** In normal adult rabbits, this antibody reacts with 30-60% of peripheral blood leukocytes. It has also been reported that the antibody stains: lymphocytes and macrophages in the spleen; lymphocytes in the thymic medulla and some cells in the cortex; Kupffer cells in the liver; and synovial lining cells in inflamed synovium. DQ? transcripts have been detected in bone marrow and spleen and, to a lesser extent, in appendix, lymph nodes, and other organs.

**细胞定位** Cell Membrane

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors