

Product datasheet

Anti-Yersinia pestis v antigen antibody [Va13] ab30808

概述

| | |
|--------------|---|
| 产品名称 | Anti-Yersinia pestis v antigen抗体[Va13] |
| 描述 | 小鼠单克隆抗体[Va13] to Yersinia pestis v antigen |
| 宿主 | Mouse |
| 特异性 | ab30808 recognises Yersinia pestis V antigen. |
| 经测试应用 | 适用于: ELISA, WB, Sandwich ELISA |
| 种属反应性 | Reacts with Yersinia pestis. Not yet tested in other species. |
| 免疫原 | Recombinant full length V-antigen of <i>Yersinia pestis</i> |
| 常规说明 | <p>Concentration of this antibody varies from lot to lot.</p> <p>Abcam is committed to meeting high standards of ethical manufacturing and as such, we will be discontinuing this product, which has been generated by the ascites method, within the next year. We are sorry for any inconvenience this may cause. If you would like help finding an alternative product, please do not hesitate to contact our scientific support team.</p> |

性能

| | |
|-------------|---|
| 形式 | Liquid |
| 存放说明 | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. |
| 存储溶液 | Preservative: 0.1% Sodium Azide Constituents: PBS, pH 7.4 |
| 纯度 | Protein G purified |
| 克隆 | 单克隆 |
| 克隆编号 | Va13 |
| 同种型 | IgG1 |

应用

Our [Abpromise guarantee](#) covers the use of **ab30808** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| 应用 | Ab评论 | 说明 |
|----------------|------|---|
| ELISA | | Use at an assay dependent dilution. |
| WB | | Use at an assay dependent dilution. Predicted molecular weight: 37 kDa. |
| Sandwich ELISA | | Use at an assay dependent dilution. Can be used as detection antibody in conjugation with recommended capture pair. |

靶标

相关性 Yersinia pestis has a number of virulence factors that enable it to survive in humans by facilitating use of host nutrients, causing damage to host cells, and subverting phagocytosis and other host defense mechanisms. Yersinia pestis is possibly involved in calcium regulation of YOP expression, which includes the export process.

细胞定位 Secreted

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors