

## Product datasheet

# Anti-Thyroxine Binding Globulin antibody [M26211] ab30832

### 概述

产品名称	Anti-Thyroxine Binding Globulin抗体[M26211]
描述	小鼠单克隆抗体[M26211] to Thyroxine Binding Globulin
宿主	Mouse
特异性	ab30832 recognises Thyroxine Binding Globulin (TBG). The clone number has been updated from (8.F.274) to (M26211) both clone numbers name the same antibody clone.
经测试应用	适用于: ELISA
种属反应性	与反应: Human
免疫原	Human Thyroxine Binding Globulin

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Add glycerol to a final volume of 50% for extra stability and aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.05% Sodium Azide Constituents: PBS, pH 7.2
纯度	Protein A purified
克隆	单克隆
克隆编号	M26211
同种型	IgG2a

### 应用

Our [Abpromise guarantee](#) covers the use of **ab30832** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
ELISA		Use at an assay dependent dilution.

### 靶标

---

---

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors