

## Product datasheet

# Anti-Sex Hormone Binding Globulin/SHBG antibody [8.F.256] ab31401

### 概述

---

产品名称	Anti-Sex Hormone Binding Globulin/SHBG抗体[8.F.256]
描述	小鼠单克隆抗体[8.F.256] to Sex Hormone Binding Globulin/SHBG
宿主	Mouse
经测试应用	适用于: ELISA, IHC-Fr, WB
种属反应性	与反应: Human
免疫原	Full length human Sex Hormone Binding Globulin protein.
常规说明	Suitable for use as the detection or conjugate (label) antibody.

### 性能

---

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Add glycerol to a final volume of 50% for extra stability and aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Constituent: PBS
纯度	Protein A purified
Primary antibody说明	Suitable for use as the detection or conjugate (label) antibody.
克隆	单克隆
克隆编号	8.F.256
同种型	IgG2a
轻链类型	kappa

### 应用

---

Our [Abpromise guarantee](#) covers the use of **ab31401** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
ELISA		Use at an assay dependent concentration.
IHC-Fr		Use at an assay dependent concentration.
WB		Use at an assay dependent concentration. Predicted molecular weight: 44 kDa.

## 靶标

<b>功能</b>	Functions as an androgen transport protein, but may also be involved in receptor mediated processes. Each dimer binds one molecule of steroid. Specific for 5-alpha-dihydrotestosterone, testosterone, and 17-beta-estradiol. Regulates the plasma metabolic clearance rate of steroid hormones by controlling their plasma concentration.
<b>组织特异性</b>	Isoform 1 and isoform 2 are present in liver and testis.
<b>序列相似性</b>	Contains 2 laminin G-like domains.
<b>细胞定位</b>	Secreted. In testis, it is synthesized by the Sertoli cells, secreted into the lumen of the seminiferous tubule and transported to the epididymis.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors