

## Product datasheet

# Anti-AU1 tag antibody (HRP) ab3402

### 概述

产品名称	Anti-AU1 tag抗体(HRP)
描述	山羊多克隆抗体to AU1 tag (HRP)
宿主	Goat
偶联物	HRP
经测试应用	适用于: ELISA, ICC, WB
免疫原	Full length native AU1 (purified): DTYRYI conjugated to KLH.

[Run BLAST with](#)

[Run BLAST with](#)

常规说明 Molar enzyme/antibody protein ratio is 4:1.

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 7.00 Constituents: 1.19% HEPES, 0.58% Sodium chloride, 0.2% BSA
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

### 应用

Our [Abpromise guarantee](#) covers the use of **ab3402** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
ELISA		1/10000 - 1/100000.
ICC		1/200 - 1/500.

应用	Ab评论	说明
WB		1/10000 - 1/30000.

## 靶标

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors