

Product datasheet

Anti-Alkaline Phosphatase antibody ab34412

概述

产品名称	Anti-Alkaline Phosphatase抗体
描述	山羊多克隆抗体to Alkaline Phosphatase
宿主	Goat
特异性	ab34412 recognises Alkaline Phosphatase (AP).
经测试应用	适用于: , ELISA
种属反应性	与反应: Cow
免疫原	Purified Alkaline Phosphatase (Calf)

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Constituent: Whole serum
纯度	Whole antiserum
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab34412** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
----	------	----

Phosphatase Activity

ELISA

应用说明	<p>ELISA: Use at an assay dependent dilution.</p> <p>APAAP (Complexes of anti-alkaline phosphatase antibody with alkaline phosphatase): Use at an assay dependent dilution. Prepare immediately prior to use and do not freeze APAAP complex.</p> <p>Not yet tested in other applications.</p>
------	---

Optimal dilutions/concentrations should be determined by the end user.

靶标

序列相似性 Belongs to the alkaline phosphatase family.

细胞定位 Cell membrane.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors