

Product datasheet

Anti-CD30-L antibody [RM153] (FITC) ab36431

1 图像

概述

产品名称	Anti-CD30-L 抗体[RM153] (FITC)
描述	大鼠单克隆抗体[RM153] to CD30-L (FITC)
宿主	Rat
偶联物	FITC. Ex: 493nm, Em: 528nm
经测试应用	适用于: Flow Cyt
种属反应性	与反应: Mouse
免疫原	Tissue/ cell preparation: CD153/CHO cells (Mouse).
常规说明	Clone RM153 is reported to block the functional binding of CD153 and can be used to inhibit the binding of murine CD30 to CD153.

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: 0.09% Sodium Azide Constituents: 1% BSA, PBS, pH 7.4
纯度	Protein G purified
Primary antibody说明	Clone RM153 is reported to block the functional binding of CD153 and can be used to inhibit the binding of murine CD30 to CD153.
克隆	单克隆
克隆编号	RM153
骨髓瘤	P3U1
同种型	IgG2b

应用

Our [Abpromise guarantee](#) covers the use of **ab36431** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
----	------	----


Flow Cyt 1/10 - 1/25.
Use 10µl of this concentration to label 10⁶ cells in 100µl.

[ab136125](#) - Rat monoclonal IgG2b, is suitable for use as an isotype control with this antibody.

靶标

功能 Cytokine that binds to TNFRSF8/CD30. Induces proliferation of T-cells.
序列相似性 Belongs to the tumor necrosis factor family.
细胞定位 Membrane.

图片



ab36431 staining CD153 in activated mouse spleen cells by Flow cytometric analysis.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors