

## Product datasheet

# Anti-Procathepsin K antibody [4b9] ab39246

### 概述

<b>产品名称</b>	Anti-Procathepsin K抗体[4b9]
<b>描述</b>	小鼠单克隆抗体[4b9] to Procathepsin K
<b>宿主</b>	Mouse
<b>特异性</b>	This antibody is specific for Procathepsin K.
<b>经测试应用</b>	<b>适用于:</b> WB, ELISA, IHC-P
<b>种属反应性</b>	<b>与反应:</b> Human
<b>免疫原</b>	N terminal His-tagged recombinant fusion protein (Human), produced in <i>E. coli</i> (36.7 kDa). It comprised 314 amino acid residues of the Human Procathepsin K and 12 additional amino acid residues (MRGSHHHHHHGS) - His Tag.

### 性能

<b>形式</b>	Liquid
<b>存放说明</b>	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
<b>存储溶液</b>	Preservative: None Constituents: 0.1M Sodium chloride, 0.05M PBS, pH 7.2
<b>纯度</b>	Protein G purified
<b>纯化说明</b>	Affinity chromatography on a column with immobilized protein G (sterile filtered).
<b>克隆</b>	单克隆
<b>克隆编号</b>	4b9
<b>同种型</b>	IgG2b

### 应用

Our [Abpromise guarantee](#) covers the use of **ab39246** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use at an assay dependent concentration. Predicted molecular weight: 27 kDa.

应用	Ab评论	说明
ELISA		Use at an assay dependent concentration.
IHC-P		Use at an assay dependent concentration.

## 靶标

<b>功能</b>	Closely involved in osteoclastic bone resorption and may participate partially in the disorder of bone remodeling. Displays potent endoprotease activity against fibrinogen at acid pH. May play an important role in extracellular matrix degradation.
<b>组织特异性</b>	Predominantly expressed in osteoclasts (bones).
<b>疾病相关</b>	Defects in CTSK are the cause of pycnodysostosis (PKND) [MIM:265800]. PKND is an autosomal recessive osteochondrodysplasia characterized by osteosclerosis and short stature.
<b>序列相似性</b>	Belongs to the peptidase C1 family.
<b>细胞定位</b>	Lysosome.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors