

Product datasheet

Penis (human) normal tissue slides ab4356

概述

产品名称	Penis (human) normal组织切片
常规说明	After excision, tissue was fixed in paraformaldehyde and processed to paraffin wax. Tissue sections were cut to a thickness of 5µm and mounted on positively charged glass slides. One tissue from each lot was stained by haematoxylin and eosin (H+E) and then microscopically examined to ensure quality. All tissue is examined / diagnosed by a licensed pathologist and is ethically obtained. Age of the donor is lot dependent. Please enquire.
经测试应用	<p>This product is an excellent control for use in immunohistochemistry (IHC) and in-situ hybridization.</p> <p>适用于: IHC-P, In situ hybridization</p>

性能

存放说明	Store at +4°C.
Slide说明	This product is an excellent control for use in immunohistochemistry (IHC) and in-situ hybridization.

应用

Our [Abpromise guarantee](#) covers the use of **ab4356** in the following tested applications. The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
IHC-P		Use at an assay dependent concentration.
In situ hybridization		Use at an assay dependent concentration.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors