


Product datasheet

Anti-GPCR GPR44 antibody ab46931

1 图像

概述

产品名称	Anti-GPCR GPR44抗体
描述	兔多克隆抗体to GPCR GPR44
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Mouse, Rat, Human 预测可用于: Chimpanzee, Macaque monkey 
免疫原	Synthetic peptide: DDSELGGAGSSRRRRTSSTA , corresponding to C terminal amino acids 329-348 of Human GPCR GPR44. Run BLAST with Run BLAST with
阳性对照	GPCR GPR44-transfected Jurkat cells.

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.02% Sodium azide Constituent: PBS
纯度	Protein G purified
克隆	多克隆
同种型	IgG

应用


Our [Abpromise guarantee](#) covers the use of **ab46931** in the following tested applications.
The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/500 - 1/1000. Detects a band of approximately 66 kDa (predicted molecular weight: 43 kDa).

靶标

功能	Orphan receptor.
序列相似性	Belongs to the G-protein coupled receptor 1 family.
细胞定位	Cell membrane.

图片



Anti-GPCR GPR44 antibody (ab46931) at 1/500 dilution + GPR44 transfected Jurkat cells

Predicted band size: 43 kDa

Observed band size: 66 kDa

[why is the actual band size different from the predicted?](#)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors