

Product datasheet

Anti-ACCN5 antibody ab48728

1 图像

概述

产品名称	Anti-ACCN5抗体
描述	兔多克隆抗体to ACCN5
宿主	Rabbit
特异性	This antibody recognizes ACCN5.
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	A region within synthetic peptide: LLDCEFFGKP CSPKDFAHVF TEYGNCFTFN HGETLQAKRK VSVSGRGLSL, corresponding to internal sequence amino acids 181-230 of Human ACCN5 Run BLAST with ExPASy Run BLAST with NCBI
阳性对照	Jurkat cell lysate

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: 0.09% Sodium azide Constituents: 2% Sucrose, PBS
纯度	Protein A purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab48728** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab 评论	说明
----	-------	----

WB Use a concentration of 1.25 µg/ml. Detects a band of approximately 48 kDa (predicted molecular weight: 57 kDa). Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

靶标

图片



Western blot - Anti-ACCN5 antibody (ab48728)

Anti-ACCN5 antibody (ab48728) at 1.25 µg/ml + Jurkat cell lysate (approximately 10µg) with 5% skimmed milk in PBS

Secondary

HRP conjugated anti-Rabbit IgG at 50,000

Predicted band size: 57 kDa

Observed band size: 48 kDa

[why is the actual band size different from the predicted?](#)

Gel concentration:12%

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors