

Product datasheet

Anti-ZNF365 antibody ab49243

2 图像

概述

产品名称	Anti-ZNF365抗体
描述	兔多克隆抗体to ZNF365
宿主	Rabbit
经测试应用	适用于: WB, ELISA, IHC-P
种属反应性	与反应: Human 预测可用于: Mouse, Rat, Rabbit, Horse, Guinea pig, Cow, Cat, Dog 
免疫原	A region within synthetic peptide: RFRSLSSLRA HLEFSHSYEE RTLLTKCSLF PSLKDTDLVT SSELLKPGKL, corresponding to N terminal amino acids 36-85 of Human ZNF365 Run BLAST with ExPASy Run BLAST with NCBI
阳性对照	HepG2 cell lysate.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: None Constituents: 2% Sucrose, PBS
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab49243** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab 评论	说明


应用	Ab 评论	说明
WB		Use a concentration of 1.25 µg/ml. Detects a band of approximately 38 kDa (predicted molecular weight: 38 kDa). Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.
ELISA		Use at an assay dependent concentration.
IHC-P		Use a concentration of 4 - 8 µg/ml.

靶标

相关性 Zinc Finger Protein 365 is a new candidate transcription factor.

细胞定位 Cytoplasm; cytoskeleton; centrosome.

图片



Anti-ZNF365 antibody (ab49243) at 1.25 µg/ml + HepG2 cell lysate at 10 µg with 5% skim milk in PBS.


Secondary

HRP conjugated anti-Rabbit IgG at 1/50000 dilution

Predicted band size: 38 kDa

Observed band size: 41 kDa

Gel concentration: 12%



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) analysis of human liver tissue labelling ZNF365 at 4-8µg/ml. Arrows indicate positively labelled hepatocytes. Magnification: 400X.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-ZNF365 antibody (ab49243)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors