

Product datasheet

Mouse monoclonal [RG-16] Anti-Rabbit IgG+IgM+IgA heavy chain (HRP) ab49626

2 References

概述

产品名称	小鼠单克隆抗体[RG-16]抗兔IgG+IgM+IgA heavy chain (HRP)
宿主	Mouse
靶标种属	Rabbit
经测试应用	适用于: ELISA, IHC-P, Dot blot
最小交叉反应	Cat, Chicken, Cow, Dog, Goat, Guinea pig, Horse, Human, Pig, Rat, Sheep, Turkey more details
免疫原	Rabbit Immunoglobulins (IgG, IgA, IgM)
表位	Recognizes an epitope located on the heavy chain of rabbit IgG, IgM, and IgA. Reduction of rabbit immunoglobulins appears to destroy the epitope.
阳性对照	Human tonsils.
偶联物	HRP

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 7.40 Preservative: 0.01% Thimerosal (merthiolate) Constituents: 0.0268% PBS, 1% BSA
纯度	IgG fraction
共轭说明	The immunoglobulin fraction of the ascites fluid is conjugated to horseradish peroxidase by protein cross-linking with 0.2% glutaraldehyde.
克隆	单克隆
克隆编号	RG-16
同种型	IgG1

应用

Our [Abpromise guarantee](#) covers the use of **ab49626** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
ELISA		
IHC-P		
Dot blot		

应用说明

Dot: 1/30000.
ELISA: 1/30000.
IHC-P: 1/100 using human tonsil tissue.

Not yet tested in other applications.
Optimal dilutions/concentrations should be determined by the end user.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors