


Product datasheet

Anti-HOXC11 antibody [HOX5J232] ab51238

1 图像

概述

产品名称	Anti-HOXC11抗体[HOX5J232]
描述	小鼠单克隆抗体[HOX5J232] to HOXC11
宿主	Mouse
经测试应用	适用于: Dot blot, WB
种属反应性	与反应: Recombinant fragment 预测可用于: Human 
免疫原	Recombinant Human HOXC11

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.05% Sodium Azide Constituents: 1% BSA, PBS, pH 7.4
纯度	Protein G purified
纯化说明	ab51238 was purified from culture supernatant of the hybridoma, cultured in a medium containing bovine IgG-depleted (approximately 95%) fetal bovine serum.
克隆	单克隆
克隆编号	HOX5J232
同种型	IgG1

应用

Our [Abpromise guarantee](#) covers the use of **ab51238** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
Dot blot		Use at an assay dependent dilution.

应用	Ab评论	说明
WB		Use at an assay dependent dilution. Predicted molecular weight: 34 kDa.


## 靶标

**功能** Sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis. Binds to a promoter element of the lactase-phlorizin hydrolase gene.

**序列相似性** Belongs to the Abd-B homeobox family.  
Contains 1 homeobox DNA-binding domain.

**细胞定位** Nucleus.

## 图片



Anti-HOXC11 antibody [HOX5J232]  
(ab51238)

**Predicted band size:** 34 kDa

**Observed band size:** 42 kDa

Western blot analysis of immunised recombinant protein.

This antibody has only been tested in WB against the recombinant fragment used as immunogen. We have no data on the detection of endogenous protein.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors