

Product datasheet

Anti-Folate Receptor 4 antibody [12A5] (FITC) ab53683

1 References 1 图像

概述

<b>产品名称</b>	Anti-Folate Receptor 4抗体[12A5] (FITC)
<b>描述</b>	大鼠单克隆抗体[12A5] to Folate Receptor 4 (FITC)
<b>宿主</b>	Rat
<b>偶联物</b>	FITC. Ex: 493nm, Em: 528nm
<b>经测试应用</b>	<b>适用于:</b> WB, IP, Flow Cyt
<b>种属反应性</b>	<b>与反应:</b> Mouse
<b>免疫原</b>	Tissue/ cell preparation (Mouse) - spleen and lymph node
<b>阳性对照</b>	Mouse spleen
<b>常规说明</b>	When this antibody is used for regulatory T cells instead of Foxp3, the cells can be sorted alive to investigate functions.

性能

<b>形式</b>	Liquid
<b>存放说明</b>	Shipped at 4°C. Store at +4°C.
<b>存储溶液</b>	Preservative: 0.001% Sodium Azide Constituents: PBS
<b>纯度</b>	Protein G purified
<b>Primary antibody说明</b>	When this antibody is used for regulatory T cells instead of Foxp3, the cells can be sorted alive to investigate functions.
<b>克隆</b>	单克隆
<b>克隆编号</b>	12A5
<b>同种型</b>	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab53683** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/1000. Predicted molecular weight: 28 kDa.
IP		1/500.
Flow Cyt		1/500.

## 靶标

相关性	Folate receptor 4 is a highly expressed cell surface marker of regulatory T cells.
细胞定位	Secreted

## 图片



Detection of cell membrane mouse folate receptor 4 in spleen of mouse Cells were stained by using PE-conjugated anti-mouse CD25, APC-conjugated anti-mouse CD4 mouse, and FITC-conjugated folate receptor4 antibody ab53683

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors