

Product datasheet

Anti-ERH antibody ab54690

2 图像

概述

产品名称	Anti-ERH抗体
描述	小鼠单克隆抗体to ERH
宿主	Mouse
经测试应用	适用于: IHC-P, WB
种属反应性	与反应: Human, Recombinant fragment
免疫原	Recombinant full length protein, corresponding to amino acids 1-105 of Human ERH

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None PBS, pH 7.2
纯度	Protein G purified
克隆	单克隆
同种型	IgG1
轻链类型	kappa

应用

Our [Abpromise guarantee](#) covers the use of **ab54690** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab评论	说明
IHC-P		Use a concentration of 3 µg/ml.

应用	Ab评论	说明
WB		Use a concentration of 1 - 5 µg/ml. This antibody has only been tested in WB against the recombinant fragment used as immunogen. We have no data on the detection of endogenous protein.

靶标


功能	May have a role in the cell cycle.
组织特异性	Expressed in all tissues examined.
序列相似性	Belongs to the E(R) family.

图片



ERH antibody (ab54690) used in immunohistochemistry at 3ug/ml on formalin fixed and paraffin embedded human lung.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-ERH antibody (ab54690)



Western blot against tagged recombinant protein immunogen using ab54690 ERH antibody at 1ug/ml. Predicted band size of immunogen is 37 kDa

Western blot - Anti-ERH antibody (ab54690)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet

- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors