

Product datasheet

Anti-ART3 antibody ab54805

2 图像

概述

产品名称	Anti-ART3抗体
描述	小鼠单克隆抗体to ART3
宿主	Mouse
经测试应用	适用于: IHC-P, WB
种属反应性	与反应: Human, Recombinant fragment
免疫原	Recombinant fragment, corresponding to amino acids 29-139 of Human ART3

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None PBS, pH 7.2
纯度	Protein G purified
克隆	单克隆
同种型	IgG2b
轻链类型	kappa

应用

Our [Abpromise guarantee](#) covers the use of **ab54805** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
----	------	----

IHC-P

WB

应用说明	WB: Use at a concentration of 1-5 µg/ml. IHC-P: Use at a concentration of 2 µg/ml.
------	---

This antibody has only been tested in WB against the recombinant fragment used as immunogen. We have no data on the detection of endogenous protein.


Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

靶标


组织特异性	Testis specific.
序列相似性	Belongs to the Arg-specific ADP-ribosyltransferase family.
细胞定位	Cell membrane.

图片



ART3 antibody (ab54805) used in immunohistochemistry at 2ug/ml on formalin fixed and paraffin embedded human testis.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-ART3 antibody (ab54805)



Western blot against tagged recombinant protein immunogen using ab54805 ART3 antibody at 1ug/ml. Predicted band size of immunogen is 38 kDa

Western blot - Anti-ART3 antibody (ab54805)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery

- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors