

## Product datasheet

# Anti-Parvin gamma antibody ab56367

### 1 图像

#### 概述

产品名称	Anti-Parvin gamma抗体
描述	小鼠单克隆抗体to Parvin gamma
宿主	Mouse
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Recombinant full length protein, corresponding to amino acids 1-332 of Human Parvin gamma

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None PBS, pH 7.2
纯度	Protein G purified
克隆	单克隆
同种型	IgG1
轻链类型	kappa

#### 应用

Our [Abpromise guarantee](#) covers the use of **ab56367** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use a concentration of 1 - 5 µg/ml. Predicted molecular weight: 37 kDa.

#### 靶标

<b>功能</b>	Probably plays a role in the regulation of cell adhesion and cytoskeleton organization.
<b>组织特异性</b>	Expressed predominantly in lymphoid organs, including spleen, thymus, lymph node, bone marrow and peripheral blood leukocytes and moderately in the digestive tract, including stomach, duodenum, jejunum, ileum, ileocecum and appendix, as well as in lung and liver. Also expressed in tumors, but at a lower level than in the corresponding normal tissues.
<b>序列相似性</b>	Belongs to the parvin family. Contains 2 CH (calponin-homology) domains.
<b>细胞定位</b>	Cell junction > focal adhesion. Cell membrane. Cytoplasm > cytoskeleton. Constituent of focal adhesions.

## 图片



Parvin gamma antibody (ab56367) at  
1ug/lane + HL-60 cell lysate at 25ug/lane.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors